

POSTER PRESENTATION

Open Access

An efficient way to analyze function of immune cells

Zhong Wang¹, Namyong Kim^{2*}, Melvin Lye²

From Society for Immunotherapy of Cancer 28th Annual Meeting
National Harbor, MD, USA. 8-10 November 2013

Cellular assays for suspension cells provide investigators the opportunity to interrogate the biology of cellular function of immune cells or immunotherapeutic drugs. Despite significant technical advances made during the last few years, working with suspension cells, such as primary blood cells, PBMCs (peripheral blood mononuclear cells), leukemia, myeloma, or lymphoma cells, cellular assays still suffer from limitations. For example, even for HEK 293 cells, a commonly used semi-adherent cell line, significant cell loss can occur during assay process. This loss can be further exacerbated when multi-step staining procedures are required.

We applied DropArray plates, which employ surface tension to retain cells on the plate surface, for our suspension cellular assays. For a variety of suspension cells, including PBMCs, AML cells, or B cells, efficient retention (50% to 80%) is achieved during multiple washes. Reproducibility and quality of the experimental results are improved with DropArray plates. Furthermore, the specific cytotoxic action of bispecific antibody on tumor cells can be also observed in real time.

Authors' details

¹Sun Yat-Sen University, Guangzhou, China. ²Curiox Biosystems, San Carlos, CA, USA.

Published: 7 November 2013

doi:10.1186/2051-1426-1-S1-P253

Cite this article as: Wang et al.: An efficient way to analyze function of immune cells. *Journal for ImmunoTherapy of Cancer* 2013 **1**(Suppl 1):P253.

Submit your next manuscript to BioMed Central and take full advantage of:

- Convenient online submission
- Thorough peer review
- No space constraints or color figure charges
- Immediate publication on acceptance
- Inclusion in PubMed, CAS, Scopus and Google Scholar
- Research which is freely available for redistribution

Submit your manuscript at
www.biomedcentral.com/submit



²Curiox Biosystems, San Carlos, CA, USA

Full list of author information is available at the end of the article